

Product Data Sheet - Evolution Passive Infrared Detector v1



Rotating Passive Infrared Detector - Battery Operated/Mains Powered

The Evo Passive Infrared detector is designed to meet the needs and legal requirements of the Alarm & Nursecall markets. Provides a high quality featured product at an economical cost.

The installation of this Passive infrared device is very simple, requiring the use of either the quick wall fix command strip system (Can be removed later with no damage to decoration) or the conventional wall bracket with two screws to fix each unit in its required position. With either of these two fixing options the device has a quick release and can be moved from room to room as required. Simple switched on/off operation with a timed delay to allow staff to exit protected area before product arms itself. Connection to all types of systems by use of a connection lead. Supplied separately for your type system when ordering.

This product is now taking the place of the every day pressure mat commonly used in the nursing home sector. Due to its location, eliminates it being damaged in the same manner as a pressure mat, and covers the whole area.

Nursecall Operation

The Evo Passive Infrared detector can be used in multiple positions to allow the detecting of the following situations - The exiting of a room. Side of bed or chair to detect client exiting. Head of bed to detect client sitting up or above chair to detect client about to stand up. The unique design projects a thin curtain beam from floor to ceiling or above the head when mounted above the bed or chair.

Features

- Quick installation, up and running in minutes
- No programming required works direct from box
- Rotating wall bracket. Point in any direction
- Thin curtain beam covers floor to ceiling
- Freedom to move to required locations
- Battery or mains powered
- Low battery warning indicator

Product code:- PIR-FS-BAT-STD

Check out this product on our webpage - http://mcd4u.com/pir_fs_bat_std.html

For more information or a quote on this versatile product please contact our sales team or follow the link above.

This product may not be suitable for use with "high risk for falls" individuals. This monitor is not a device used for falls prevention - Its is used to detect a potential fall and cannot prevent a fall.